

Detecting GMOS via PCR & Gel Electrophoresis

By Darren Zheng

11th Grade

Introduction

- Genetic Engineering
 - Encode beneficial traits into plants and animals
- 448 million acres of GM crops globally in 2014
 - 82% of Soy
 - 68% of Cotton
 - 30% of Corn
 - 25% of Canola

Background Info

- Genetically Modified Organism
 - “Unnatural” changes
- Transgenes – genes from another organism
- 2 Important Example
 - *Bacillus thuringiensis* – Toxic to Corn Borers
 - Glyphosate-resistance gene – protects from common herbicides
 - Higher yield, lower environmental impact, etc.

Objective

- To determine if certain “Non-GMO verified” claims are valid (35S Promotor and 225 NOS)
 - Triscuit Thin Crisps
 - Soybeans (Generic Brand)
 - Banza Spagetti
 - Popcorn Cakes
 - Peanuts (Generic Brand)
 - Post Great Grains Cereal



Foods



Hypothesis

- I believe that the non-GMO verified claims on the foods that have them are true.
 - Foods w/ “Non-GMO”
 - Cereal, Triscuit, Spagetti, Corn
 - Foods w/o Labels
 - Peanuts & Soybeans

Procedure

- Isolation
 - Extract samples of DNA from food products
- Amplification
 - Use PCR to replicate the DNA for larger pool
- Gel Electrophoresis
 - Look for the specific DNA bands to determine if GMO

Materials

- Mortar & Pestle
- Non-GMO Grain Sample (The Control)
- Pipettes & Distilled Water
- Water Bath & Assorted Tubes
- Insta-Gene Matrix (Microscopic Beads)
- Scale
- Centrifuge

Isolation

1. Pulverize the food sample in a sterile mortar & pestle into a fine powder
 - A. Weigh out ~1 gram and mix with 5 ml of water
2. Pipet 50 μ l of slurry into a tube with 500 μ l of Insta-Gene
3. Repeat previous steps for each food sample & control
4. Water Bath @ 95° C for 5 min, then centrifuge

Materials

- PCR Thermal Cycler
- Plant PSII (green) & GMO (red) primers
- Master Mix
- GMO-positive sample (Another Control)
- Pipettes & Distilled Water
- Ice-bath & Assorted Tubes

Amplification

1. Pulse-spin all tubes & template (50 μ l); Avoid Supernatant
2. Dilute primers with the Master Mix
 - A. Ratio of 11 μ l primer : 550 μ l master mix; Store on Ice
3. Create the following mixtures via pipette in PCR tubes:
 - A. 20 μ l plant + 20 μ l food sample
 - B. 20 μ l GMO + 20 μ l food sample
4. Repeat for all test foods and controls; Made 2 of each

Amplification (pt. 2)

5. Place the PCR tubes in the Thermal Cycler



A. 94° C @ 2 min

B. Begin 40 Cycles:

A. 94° C @ 1 min

B. 59° C @ 1 min

C. 72° C @ 2 min

C. 72° C for 10; Hold @ 4° C

Materials

- Agarose Gel (Made from Powder)
- 1X TAE buffer (Tris-acetate-EDTA) & Methylene Blue Dye
- Orange Loading Dye & Molecular Weight Ruler
- Pipettes & Distilled Water & Assorted Tubes/Flasks
- Transilluminator (LightBox)
- Gel Electrophoresis Chamber & Equipment
- Weighting Dishes/Ziplock Bags

Gel Electrophoresis

1. Prepare 3% agarose gels with 6 wells each & place in station
 - A. Set up gel electrophoresis station & pour TAE buffer in
2. Remove PCR tubes and place in a tray
3. Pipette 10 μ l of orange loading dye into each tube
 - A. Add 40 μ l to 1000 μ l of MWR & mix well
4. Load 20 μ l of each sample into individual lanes
 - A. Requires 1 MWR in each gel




Pictures



Gel Electrophoresis (pt. 2)

5. Run the gels until ~80% orange front (150 V @ 45 min)
6. Stain with Methylene blue dye overnight in weighing boats
 - A. 1X concentration (Diluted from 500X)/Ziplock bags
7. Wash off gel and soak in distilled water
 - A. Be careful of dehydration
8. Repeat until bands are visible under light box
9. Record the data and analyze

Pictures




Diagram

Note: Control Refers to POSITIVE GMO

	Gel #1	Gel #2	Gel #3	Gel #4	Gel #5	Gel #6	Gel #7
Lane #1	MWR	MWR	MWR	MWR	MWR	MWR	MWR
Lane #2	Cereal GMO	Spaghetti GMO	Corn GMO	Triscuit GMO	Peanuts GMO	Soy GMO	Grain GMO
Lane #3	Cereal GMO	Spaghetti GMO	Corn GMO	Triscuit GMO	Peanuts GMO	Soy GMO	Grain GMO
Lane #4	Cereal Plant	Spaghetti Plant	Corn Plant	Triscuit Plant	Peanuts Plant	Soy Plant	Grain Plant
Lane #5	Cereal Plant	Spaghetti Plant	Corn Plant	Triscuit Plant	Peanuts Plant	Soy Plant	Grain Plant
Lane #6	Control GMO	Control GMO	Control GMO	Control GMO	Control Plant	Control Plant	Control Plant

Pictures



True/False Chart

	Gel #1	Gel #2	Gel #3	Gel #4	Gel #5	Gel #6	Gel #7
Lane #1	MWR	MWR	MWR	MWR	MWR	MWR	MWR
Lane #2	Cereal GMO	Spaghetti GMO	Corn GMO	Triscuit GMO	Peanuts GMO	Soy GMO	Grain GMO
Lane #3	Cereal GMO	Spaghetti GMO	Corn GMO	Triscuit GMO	Peanuts GMO	Soy GMO	Grain GMO
Lane #4	Cereal Plant	Spaghetti Plant	Corn Plant	Triscuit Plant	Peanuts Plant	Soy Plant	Grain Plant
Lane #5	Cereal Plant	Spaghetti Plant	Corn Plant	Triscuit Plant	Peanuts Plant	Soy Plant	Grain Plant
Lane #6	Control GMO	Control GMO	Control GMO	Control GMO	Control Plant	Control Plant	Control Plant

Conclusion

- My hypothesis was partially correct.
- Determined to be non-GMO
 - Cereal, Spaghetti, Corn, Triscuit, Peanuts
- Soy discarded due to Error, etc.
- The controls worked as expected for the majority
- Plant Primer?

Limitations/Errors

- Voltage control
- Sensitivity of agarose gels during formation
- Staining issues with band variation
- Lack of time/resources to process more trials
- Potential Contamination in Mortar/Pestle
- Inconsistency between separate gel runs

Further Inquiry

- More trials and larger sample sizes
- More types of food
- What is the threshold for non-GMO certified?
- Calculate Retardation factor to normalize results
- Perform horizontal gene transfer in plant growth
- CRISPR research

Works Cited

“Genetically Modified Organism.” *Wikipedia*, Wikimedia Foundation, 23 Feb. 2019, en.wikipedia.org/wiki/Genetically_modified_organism.

“Horizontal Gene Transfer.” *Welcome -Antimicrobial Resistance Learning Site For Veterinary Students*, amrls.cvm.msu.edu/microbiology/molecular-basis-for-antimicrobial-resistance/acquired-resistance/acquisition-of-antimicrobial-resistance-via-horizontal-gene-transfer.

Valle. “Detecting Genetically Modified Food by PCR.” *Carolina Biological Supply: World-Class Support for Science & Math*, 3 Nov. 2015, www.carolina.com/pcr-kits/detecting-genetically-modified-food-by-pcr/FAM_211366.pr.

Questions



Extra